

Main

May 11, 2025

Connected to venv (Python 3.12.4)

```
[ ]: import utils.DataLoader as dl
import utils.BiopackAnalysis as ba
import heartpy as hp
import matplotlib.pyplot as plt
import numpy as np

path = "C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony/ai.acq"

#badanie = dl.load_data(path)
#ba.draw_plots(badanie)

badania = dl.mass_load("C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony")

ba.draw_plots(badania[3])

data = badania[0]
sample_rate = 1000 # Hz - dostosuj do swojego sygnału

_, _ = ba.heartrate_analysis(data, sample_rate)
badania_filtr = data.copy()
badania_filtr['PPG (Pulse)'] = ba.bandpass_ppg(badania_filtr['PPG (Pulse)'],
↪sample_rate)
badania_filtr['EMG (5 - 500 Hz)'] = ba.bandpass_emg(badania_filtr['EMG (5 - 500_
↪Hz)'], sample_rate)
badania_filtr['EMG (5 - 500 Hz) 1'] = ba.bandpass_emg(badania_filtr['EMG (5 -_
↪500 Hz) 1'], sample_rate)
badania_filtr['Respiration'] = ba.bandpass_resp(badania_filtr['Respiration'],
↪sample_rate)
ba.draw_plots(badania_filtr)
_, _ = ba.heartrate_analysis(badania_filtr, sample_rate)

_ = ba.respiratory_analysis(badania_filtr, sample_rate)
```

Znalezienie pliki .acq:

1. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\ai.acq
2. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\hak.acq

3. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\klasik.acq
4. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\lif.acq
5. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\nothing.acq
6. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\reference.acq
7. C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\small.acq

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura -
Projekt/legia/jhony\ai.acq

C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\ai.acq

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: PPG (Pulse), Jednostka: mV

Nazwa kanału: Respiration, Jednostka: mV

	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	-0.013733	0.002747	0.054810	0.028381
0.001	-0.026550	0.006714	0.053955	0.028076
0.002	-0.029602	0.008545	0.053101	0.027771
0.003	-0.025024	0.006104	0.052185	0.027466
0.004	-0.021973	0.004272	0.051270	0.027161

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura -
Projekt/legia/jhony\hak.acq

C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\hak.acq

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: PPG (Pulse), Jednostka: mV

Nazwa kanału: Respiration, Jednostka: mV

	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	-0.044556	-0.002136	0.099670	0.192566
0.001	-0.033875	-0.000916	0.099060	0.191956
0.002	-0.018616	-0.003357	0.098450	0.191345
0.003	-0.006714	0.000916	0.097900	0.190735
0.004	-0.002136	0.002747	0.097351	0.190125

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura -
Projekt/legia/jhony\klasik.acq

C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\klasik.acq

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV

Nazwa kanału: PPG (Pulse), Jednostka: mV

Nazwa kanału: Respiration, Jednostka: mV

	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	0.011597	-0.009460	-0.038940	-0.088501
0.001	0.008240	-0.004578	-0.040466	-0.087891
0.002	0.007019	0.000305	-0.041931	-0.087280

0.003	0.005188	0.006409	-0.043457	-0.086975
0.004	0.003662	0.004578	-0.044983	-0.086365

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\lif.acq
C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\lif.acq
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: PPG (Pulse), Jednostka: mV
Nazwa kanału: Respiration, Jednostka: mV

	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	0.013428	-0.001221	0.102844	0.332336
0.001	0.006714	-0.002747	0.101807	0.332031
0.002	0.001526	-0.007935	0.100769	0.331726
0.003	-0.004272	-0.005798	0.099854	0.331421
0.004	-0.002747	0.000305	0.098999	0.330811



Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\nothing.acq
C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\nothing.acq
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: PPG (Pulse), Jednostka: mV
Nazwa kanału: Respiration, Jednostka: mV



	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	-0.005798	-0.005798	-0.152283	-0.094604
0.001	-0.005798	-0.004272	-0.152588	-0.094604
0.002	-0.009460	-0.002136	-0.152893	-0.094604
0.003	-0.010071	-0.003357	-0.153198	-0.094604
0.004	-0.004883	-0.003967	-0.153503	-0.094604

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\reference.acq
C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\reference.acq
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
Nazwa kanału: PPG (Pulse), Jednostka: mV
Nazwa kanału: Respiration, Jednostka: mV

	EMG (5 - 500 Hz)	EMG (5 - 500 Hz) 1	PPG (Pulse)	Respiration
Czas [s]				
0.000	-0.034485	-0.013733	-0.001831	0.073547
0.001	-0.047302	-0.012512	-0.002686	0.073242
0.002	-0.046387	-0.013733	-0.003540	0.072937
0.003	-0.041809	-0.012817	-0.004456	0.072327
0.004	-0.027466	-0.019836	-0.005310	0.072021

Wczytywanie danych z pliku: C:/Users/jan/Documents/Aparatura -
 Projekt/legia/jhony\small.acq
 C:/Users/jan/Documents/Aparatura - Projekt/legia/jhony\small.acq
 Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
 Nazwa kanału: EMG (5 - 500 Hz), Jednostka: mV
 Nazwa kanału: PPG (Pulse), Jednostka: mV
 Nazwa kanału: Respiration, Jednostka: mV
 EMG (5 - 500 Hz) EMG (5 - 500 Hz) 1 PPG (Pulse) Respiration
 Czas [s]
 0.000 -0.003357 -0.003052 0.100220 -0.098877
 0.001 0.000916 -0.004578 0.100037 -0.098877
 0.002 0.023804 -0.002747 0.099915 -0.099182
 0.003 0.030212 -0.002136 0.099731 -0.099487
 0.004 0.029907 -0.002441 0.099609 -0.099792





```
c:\Users\jan\Documents\Aparatura - Projekt\venv\Lib\site-
packages\heartpy\analysis.py:522: UserWarning: Short signal.
```

```
-----Warning:-----
```

```
too few peak-peak intervals for (reliable) frequency domain measure computation,
frequency output measures are still computed but treat them with caution!
```

HF is usually computed over a minimum of 1 minute of good signal. LF is usually computed over a minimum of 2 minutes of good signal. VLF is usually computed over a minimum of 5 minutes of good signal. The LF/HF ratio is usually computed over minimum 24 hours, although an absolute minimum of 5 min has also been suggested.

For more info see:


Shaffer, F., Ginsberg, J.P. (2017), An Overview of Heart Rate Variability Metrics and Norms.



Task Force of Pacing and Electrophysiology (1996), Heart Rate Variability, in: European Heart Journal, vol.17, issue 3, pp354-381




This warning will not repeat

```
warnings.warn(msg, UserWarning)
```

```
c:\Users\jan\Documents\Aparatura - Projekt\venv\Lib\site-  
packages\heartpy\visualizeutils.py:119: UserWarning: FigureCanvasAgg is non-  
interactive, and thus cannot be shown  
    fig.show()  
c:\Users\jan\Documents\Aparatura - Projekt\venv\Lib\site-  
packages\heartpy\visualizeutils.py:299: UserWarning: FigureCanvasAgg is non-  
interactive, and thus cannot be shown  
    fig.show()  
  
bpm: 76.73  
ibi: 781.99  
sdnn: 34.04  
sdsd: 15.80  
rmssd: 27.05  
pnn20: 0.47  
pnn50: 0.05  
hr_mad: 17.00  
sd1: 19.12  
sd2: 44.28  
s: 2660.24  
sd1/sd2: 0.43  
breathingrate: 0.13  
vlf: 55.54  
lf: 488.44  
hf: 145.23  
lf/hf: 3.36  
p_total: 689.20  
vlf_perc: 8.06  
lf_perc: 70.87  
hf_perc: 21.07  
lf_nu: 77.08  
hf_nu: 22.92
```










bpm: 76.73
ibi: 781.96
sdnn: 33.51
sdsd: 15.44
rmssd: 24.41
pnn20: 0.37
pnn50: 0.05
hr_mad: 18.50

sd1: 17.26
sd2: 44.27
s: 2400.14
sd1/sd2: 0.39
breathingrate: 0.13
vlf: 57.99
lf: 517.18
hf: 146.43
lf/hf: 3.53
p_total: 721.59
vlf_perc: 8.04
lf_perc: 71.67
hf_perc: 20.29
lf_nu: 77.93
hf_nu: 22.07





Średnia częstość oddechu: 16.82 oddechów/min

Średni odstęp między oddechami: 3.57 s

Odchylenie standardowe odstępów: 0.51 s

Współczynnik zmienności: 14.41 %

